

 || YUCERA



8 materials can be placed

# YRC-6X

5 Axis Wet Milling Machine




Chairside processing for instant repair in 2 hours

## Product Parameters

Specification	Parameter	Specification	Parameter
Equipment size	535*655*720mm	Processing methods	Five-axis linkage, wet milling
Weight	170kg	Burs specifications	Special props for handle diameter 4mm
Rotating axis operating angle	A: 360° B: +20°	Main axis speed	0-60,000rpm
Total power	1200w	Working pressure	4.5-7.5bar (no water, no gasoline )
Cutting precision	0.01mm	Tool magazine capacity	10

## Machinable products



## Machinable material



## Product Features




Five-axis linkage



10 tool magazines



Integrated water level monitoring



Compatible with various materials